FARMALABOR tech

DUROMETRO

HARDNESS TABLETS MACHINE

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Index

1.0 - Introduction	14
1.1 - Use of handbook	14
2.0 - General description	14
3.0 - Instructions and maintenance	14
3.1 - Scope of device	14
3.2 - Technical description	15
3.3 - How to use tablet hardness tester - procedure	15
3.4 - Periodical check of spring adjustment	15
4.0 - Rules applied	15
4.1 - Harmonized rules applied	15
4.2 - List of materials used	15
4.3 - Residual risks prevention (en 414)	16
4.4 - Producer's Address For Information - maintenance	16
5.0 - Technical data	17
5.1 - Size - weight	17
5.2 - Technical data about spring adjustment	17
5.3 - Mechanical test report	17
5.4 - Declaration of compliance	17
5.5 - Maintenance check sheet	18
5.6 - Device mechanical drawings: dis. 001, dis.002	19 - 20
Remarks	21 - 23

1.0 - Introduction

Thanks for buying our equipment and for trusting our company. The items described in this handbook have been designed and produced to fully satisfy the needs of our clients, with particular attention to the choice of materials and to their quality. All equipment has been strictly checked before shipment.

1.1 - Use of handbook

In this handbook you will find a comprehensive description of the capsule filling system, its technical characteristics and maintenance instructions. Keep it at hand for immediate use. If you lose it, please do not hesitate to ask another copy directly to FARMALABOR s.r.l. Before using your filling system please read this handbook carefully. The right functioning of your equipment depends on the right application of the instructions in this handbook. In addition, as the capsule filling system is a manual tool, the best of productivity will be reached after a period of training.

2.0 - General description

The Farmalabor Tech “Tablet Hardness Tester” is the essential tool to be used to determine the hardness of the tablets.

The hardness tester components are in chromed brass.

Maximum care have been taken to ensure optimal operation during normal conditions of use and to prevent injury.

3.0 - Instructions and maintenance

3.1 - Scope of device

The Farmalabor Tech device “Tablet Hardness Tester” has to be used only to determine the hardness of the tablets for pharmaceutical, herbal and food use.

The device can test tablets with diameter between 3 mm Ø and 24 mm Ø.

The graduated scale is between 0 daN and 15 daN.

The value is expressed in daN.

3.2 - How to use tablet hardness tester - procedure

- a) Put the tablet inside the front ferrule (10) (ref. to Pic. 2 - Pic. 002).
- b) Turn the control knob (1) to slightly press the tablet.
- c) Place the indicator pointer (5) on the 0 value of the graduated scale.
- d) Start tablet compressions, turning the control knob (1) until tablet breaks.
- e) Read measure performed (expressed in daN) on the graduated scale (11)

Note: Every imprints on graduated scale (11) correspond to 1 daN.
(Pic. 001 and Pic. 002)

3.3 - Periodical check of spring adjustment

You need to replace the spring (7) after a working cycle of one year, in order to ensure the most accurate reading on the graduated scale of the Device.

4.0 - Harmonized rules and informations

4.1 - Harmonized rules applied

EN 292 – 1

EN 292 – 2

EN 414

UNI EN 10204 3.1 B (FOR SPRING)

4.2 - List of materials used

1. Control knob: bakelite

2. Tapping screw (for control knob): stainless steel AISI 304

3. Rear ferrule: stainless steel AISI 304

- 4. Spring presser:** stainless steel AISI 304
- 5. Indicator:** stainless steel AISI 304
- 6. Socket screw:** stainless steel AISI 304
- 7. Spring:** stainless steel AISI 304
- 8. Body:** stainless steel AISI 304
- 9. Push rod:** stainless steel AISI 304
- 10. Front ferrule:** stainless steel AISI 304
- 11. Graduated scale:** stainless steel AISI 304
- 12. Linchpin:** stainless steel AISI 304
- 13. Socket screw:** stainless steel AISI 304
- 14. Lock worm screw:** stainless steel AISI 304

4.3 - Residual risks prevention

- When designing the “Tablet Hardness Tester” the residual risks were analysed in detail.
- The spring for operating the push rod (9) is completely protected and must be replaced after a year of service.
- All materials used, as specified in Paragraph “LIST OF MATERIALS”, are not dangerous for the operator’s health.
- All sharp edges on the components were eliminated in order to avoid any danger for the operator.
- All locking screws are made of stainless steel AISI 304.
- An accurate calibration of the device at producer’s site totally eliminates all residual risks.

4.4 - Producer’s address for information and maintenance

Farmalabor Tech a division of Farmalabor s.r.l.

Via Oberdan, 52

76012 – Canosa di Puglia – BT- Italy

Tel. +39 0883 611301

Fax 800 085 708

5.0 - Technical data

5.1 - Size and weight

Size: external diameter 40 x 250 mm

Weigh: 0,6 kg

5.2 - Technical data about spring adjustment

Farmalabor Tech, a division of Farmalabor S.r.l., certifies that the inspection tests carried out on the spring (particular N. 7 – Pic. N. 001 have reported the calibration values listed below:

Extension: 10 mm = 4,25 kg

Extension: 20 mm = 8,5 kg

Spring has been made using following material:

Steel C 98

Treatment: phosphate coating

5.3 - Mechanical test report

DEVICE: **Hardness Tablets Machine.**

The device was verified mechanically with regard to the spring operation and to calibration needed for measurement of the tablets hardness with a diameter ranging from 3 mm Ø to 24 mm Ø.

Furthermore all locking devices were rigorously verified.

5.4 - Declaration of compliance

We hereby certify that **Tablet Hardness Tester** complies with the harmonized European Standards and with the legislation in force.

All the components were constructed with non-toxic materials and they do not represent danger to operator's health.

Applied regulations and laws list:

EN 292-1

EN 292-2

EN 414

EN 349

UNI EN 10204 3.1 B

5.5 - Maintenance check sheet

Device: tablet hardness tester

Production year:

Maintenance n°:

Data:

Type of maintenance:

Parts replaced:

Reason for failure:

The device has been checked at the customer's site:

the device has been verified at the manufacturer's facility:

Released certificate n°:

Tester Signature

Farmalabor S.r.l.

Via Oberdan, 52

76012 – Canosa di Puglia – BT- Italy

Tel. +39 0883 611301

Fax. 800 085708

www.farmalabor.it


tech@farmalabor.it

Picture 001

1. Graduated scale

Picture 002

- 1. Push rod
- 2. Front ferrule
- 3. Body
- 4. Graduated scale
- 5. Spring presser
- 6. Tapping screw (for control knob)
- 7. Control knob
- 8. Rear ferrule
- 9. Lock worm screw
- 10. Indicator
- 11. Socket screw
- 12. Linchpin
- 13. Spring

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