Bagnomaria elettrico/ Water baths

**BMTECH 4.5<sup>®</sup>**

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## 0 SUMMARY

|       |   |    |
|-------|---|----|
| 1     | Introduction .....  | 4  |
| 1.1   | Use of handbook .....   | 4  |
| 1.1.1 | Symbology and terminology used in the use and maintenance manual..... | 4  |
| 1.2   | Statement and guaranty of the manufacturer .....                      | 5  |
| 1.3   | Warranty loss.....  | 5  |
| 2     | Field of application .....  | 5  |
| 3     | General description .....   | 6  |
| 3.1   | Generalities.....   | 6  |
| 3.1.1 | Technical features .....  | 7  |
| 3.1.2 | Contents of the package .....   | 7  |
| 3.2   | Identification.....   | 8  |
| 4     | Installation of the machine .....                                     | 8  |
| 4.1   | Environmental conditions .....  | 8  |
| 4.2   | Positioning and adjustment .....                                      | 8  |
| 4.3   | Preparations energy source .....                                      | 9  |
| 5     | Directions and maintenance .....                                      | 9  |
| 5.1   | Control .....   | 9  |
| 5.2   | Setting up of the device .....  | 9  |
| 5.2.1 | Electrical connections.....   | 9  |
| 5.2.2 | Trial start up .....  | 9  |
| 5.2.3 | Switching on the device.....  | 10 |
| 5.2.4 | Set phase .....   | 10 |
| 5.2.5 | Normal stop procedure .....   | 10 |
| 5.3   | Long time inactivity.....   | 10 |
| 5.4   | Right and wrong use .....   | 10 |
| 5.4.1 | Unintended uses.....  | 11 |
| 5.5   | Residual risks.....   | 11 |
| 5.6   | Shipment.....   | 11 |
| 5.7   | Maintenance and reparation .....                                      | 11 |
| 5.7.1 | Assembly/ disassembly.....  | 11 |

|       |  |    |
|-------|--|----|
| 5.7.2 | Cleaning suggestions .....                             | 11 |
| 6     | Electrical wiring diagram .....                        | 12 |
| 7     | Problem resolution.....                                | 13 |
| 8     | Scrapping and material disposal.....                   | 13 |
| 9     | List of applying legal abd technical specifiction..... | 14 |

# 1 INTRODUCTION

Thanks for buying our equipment and for trusting our company. The items described in this handbook have been designed and produced to fully satisfy the needs of our clients, with particular attention to the choice of materials and to their quality. All equipment has been strictly checked before shipment.

## 1.1 USE OF HANDBOOK

This handbook you will find a comprehensive description of the complete equipment, the technical characteristics and the instructions for use and maintenance.

**The manual must be kept in a suitable place.**

In case of loss or deterioration, the replacement documentation must be requested directly to FARMALABOR S.R.L.

**Before using your filling system please read this handbook carefully.**

The right functioning of your equipment depends on the right application of the instructions in this handbook.

### 1.1.1 SYMBOLOGY AND TERMINOLOGY USED IN THE USE AND MAINTENANCE MANUAL

|   |  |
|---|--|
|    | <p><b>Note</b></p> <p>Indicates useful information for consulting the manual and for the machine to work properly.</p>   |
|   | <p><b>Personal protective equipment</b></p> <p>Indicates the obligation to use protective gloves</p>   |
|  | <p><b>Attention</b></p> <p>Indicates instructions that only if correctly performed avoid situations of risk to people, informs about the dangers and how to avoid them, suggests behavioural procedures.</p> |
|  | <p><b>Electrical hazard</b></p> <p>Indicates the danger of electrocution</p>   |
|  | <p><b>Hot surfaces</b></p> <p>Indicates the presence of hot surfaces</p>   |
|  | <p><b>Prohibition</b></p> <p>Indicates actions and procedures that are prohibited.</p>   |

## 1.2 STATEMENT AND GUARANTY OF THE MANUFACTURER

All products are carefully checked, tested before shipment and are guaranteed by the company Farmalabor s.r.l. for a period of twelve (12) months from the date of delivery, except for parts subject to wear and tear. The warranty does not cover costs for labour and refers only to parts which are found defective for verified cause, attributed to manufacturing fault or poor quality of the material, with the exclusion of any right for the buyer to request any compensation for damages from any cause. Farmalabor s.r.l. or the authorized dealer who made the sale is bound to provide labour for warranty repair, the cost of which will be at the expense of the purchaser.

The warranty does not cover any damage occurred during transport, those caused by incorrect installation and /or improper use (made not in accordance with what is indicated in the user manual of the product), from inexperience and /or poor supervision.

The warranty does not include the replacement of the product and/or extension in case of any failure occurred.

## 1.3 WARRANTY LOSS

The company FARMALABOR S.R.L. declines all responsibility in case of:

- Use of the equipment contrary to the laws concerning Safety and Accident Prevention;
- Failure or incorrect observance of the instructions provided in the manual;
- Unauthorized modifications to the equipment;
- Use by untrained personnel;
- Shortage or insufficient maintenance;
- Use of spare parts not original or not indicated by the manufacturer;
- Customer outstanding with payments.

The intended use and the configurations foreseen by the machine are the only ones allowed by the manufacturer.

**Do not attempt to use the machine in disagreement with the instructions provided.**

## 2 FIELD OF APPLICATION

The water bath is an essential equipment for the production of creams, liquids and semi-liquids that need heating for the melting of the substances used.

### 3 GENERAL DESCRIPTION

#### 3.1 GENERALITIES

The BMTECH 4.5 water bath has a capacity of 4.5 l. Entirely made of stainless steel AISI 304, it is equipped with a thermostat with probe placed in contact with the tank to regulate the operating temperature.

The instrument must be placed in a covered area not exposed to the weather.

The room where the instrument is used must be lighted sufficiently to allow for operation without risks due to areas in shadow.


The accepted temperatures for normal functioning is of a minimum of +5°C a maximum of +40°C with maximum relative humidity of 80% for temperatures up to +31°C with linear reduction 50% at +40°C.

Fig. 1: Controls

1. Handles;
2. Controls;
3. Tank.

## 3.1.1 TECHNICAL FEATURES

|   |                             |
|---|-----------------------------|
| <i>Model</i>                                    | <b>BMTECH 4.5</b>           |
| <i>Tank capacity</i>                            | <b>4,5 l</b>                |
| <i>Temperature</i>                              | <b>From +30°C to +100°C</b> |
| <i>Tank precision with water</i>                | <b>± 1,5</b>                |
| <i>Voltage</i>                                  | <b>230V /50 Hz</b>          |
| <i>Power supply</i>                             | <b>1000 W</b>               |
| <i>Protection class</i>                         | <b>IP 46</b>                |
| <i>Internal size <math>\phi \times H</math></i> | <b>22 cm x 15 cm</b>        |
| <i>External size <math>\phi \times H</math></i> | <b>31 cm x 26 cm</b>        |
| <i>Weight</i>                                   | <b>5,2 kg</b>               |

## 3.1.2 CONTENTS OF THE PACKAGE

The package include:

- Water baths equipment;
- Power supply cable;
- Use and maintenance Handbook.

### 3.2 IDENTIFICATION

The EC marking on the lower part of the device, shows the main features, Fig. 2.

Fig. 2: EC marking

## 4 INSTALLATION OF THE MACHINE

### 4.1 ENVIRONMENTAL CONDITIONS

The machine must have sufficient natural light and must be equipped with adequate artificial light and ventilation to safeguard the safety and health of the operator.

The lighting of the room must comply with the laws in force in the country where the machine is installed, must guarantee a good visibility of the product.

### 4.2 POSITIONING AND ADJUSTMENT

Place the device on the working plan, providing sufficient space for use, maintenance and use of any peripheral equipment.

### 4.3 PREPARATIONS ENERGY SOURCE

Set up a bipolar single-phase outlet with earth ground at the work site


**We accept no responsibility for accidents resulting from bad connections.**

## 5 DIRECTIONS AND MAINTENANCE

### 5.1 CONTROL

The device is equipped with the following controls buttons, Fig. 3:

**Fig. 3: Controls**

- (A): Temperature regulation °C;
- (B): Main switch;
- (C): Active heating indicator.

### 5.2 SETTING UP OF THE DEVICE

#### 5.2.1 ELECTRICAL CONNECTIONS

The instrument has a single external electrical supply point placed on the back of the same.

#### 5.2.2 TRIAL START UP

Before regular and continuous use the device, check that it works perfectly by conducting at least one full trial cycle.

### 5.2.3 SWITCHING ON THE DEVICE

- Plug in the power plug into the socket placed on the back of the device;
- Plug in the supply line in at the mains;
- Turn to position "1" the main switch, Fig.

### 5.2.4 SET PHASE

- Turn the thermostat knob (A) clockwise to set the required temperature;
- The indicator light (C), Fig. 3, will light up and indicate the heating phase of the bath;

### 5.2.5 NORMAL STOP PROCEDURE

- Turn to position "0" the thermostat knob, Fig. 3;
- The warning light (C), Fig. 3, if turned on, goes off;
- Set the main switch (B) to position "0";
- Remove the necessary substances from the instrument;

#### **WARNING: to handling hot objects**

- Disconnect the device from the power supply line.

## 5.3 LONG TIME INACTIVITY

If you do not use the equipment for a long time, you have to keep it cleaned and protected from dust, humidity or any other external element that might damage any of its components.

#### **Disconnect the device from the power supply.**

## 5.4 RIGHT AND WRONG USE

The equipment must be used in full compliance with the instructions of this Handbook; any measurement of the machine might damage the system or put in danger persons or things.

### 5.4.1 UNINTENDED USES

The device must not be switched on with the empty container.

**Do not wash the device with jets of water.**

**FARMALABOR S.R.L. declines any responsibility resulting from any wrong use of the equipment.**

### 5.5 RESIDUAL RISKS

During operation, the user's limbs must not come into contact with the hot surfaces of the device and with the substances present in the tank.

Use gloves or tools suitable for handling hot containers.

### 5.6 SHIPMENT

All parts of the equipment have been checked before shipment and are supplied to customers ready for use. Upon delivery of the equipment, use this handbook description to check that nothing has been damaged during transport and that no part is missing.

### 5.7 MAINTENANCE AND REPAIRATION

**Before carrying out maintenance and / or repair operations, disconnect the device from the power supply, and make sure that the device is cold.**

Maintenance can be carried out in safety by expert personnel trained in preparation and maintenance of the instrument and authorised for the use.

#### 5.7.1 ASSEMBLY/ DISASSEMBLY

For repair operations, contact FARMALABOR S.R.L.

#### 5.7.2 CLEANING SUGGESTIONS


Carry out cleaning operations with the instrument switched off and unplugged from the mains.

We recommend daily cleaning and in any case at the end of each working day and / or work session.

Use a damp cloth and suitable neutral detergents.

Dry carefully, using a soft cloth.

## 6 ELECTRICAL WIRING DIAGRAM

| Denomination | Description           |
|--------------|-----------------------|
| <i>F1-F2</i> | Panel fuses           |
| <i>TS</i>    | Thermostat            |
| <i>IG</i>    | Main switch           |
| <i>HR</i>    | Electrical resistance |
| <i>HL-GR</i> | Led green switch      |
| <i>HL-RD</i> | Led red               |

## 7 PROBLEM RESOLUTION

| <b><i>Problem</i></b>               | <b><i>Cause</i></b>   | <b><i>Resolution</i></b>   |
|-------------------------------------|---|--|
| <b>The device doesn't switch on</b> | <ol style="list-style-type: none"> <li>1. The device is non property connected to the mains;</li> <li>2. There are problems with the mains line.</li> </ol>   | <ol style="list-style-type: none"> <li>a. Connect the device to the mains;</li> <li>b. Check the supply line;</li> </ol>   |
| <b>The device doesn't heat</b>      | <ol style="list-style-type: none"> <li>1. The device switched off;</li> <li>2. The set temperature is inferior to the liquid one in the tank;</li> <li>3. The thermostat is broken;</li> <li>4. The fuse is broken</li> </ol> | <ol style="list-style-type: none"> <li>a. Turn on the device;</li> <li>b. Increase the set temperature;</li> <li>c. Contact Farmalabor customer service;</li> <li>d. Contact Farmalabor customer service;</li> </ol> |

## 8 SCRAPPING AND MATERIAL DISPOSAL

The machine consists mainly of steel and metal parts, of plastics and rubber and mixed components.

Disposal must be carried out by delivering the whole machine or its parts to a specific company, specialized and authorized according to the laws in force, on the waste production site.

The machine and the individual parts (electrical and electronic in particular) comply with the RoHS directives 2002/95 / EC, RAEE 2002/96 / EC and 2003/108 / EC on the waste of electrical and electronic equipment; it excludes the presence of substances considered harmful such as cadmium, mercury, lead, hexavalent chromium, polybrominated biphenyls which can cause allergies or other diseases.

## 9 LIST OF APPLYING LEGAL AND TECHNICAL SPECIFICATION

The equipment has been conceived and manufactured, conforming to the following legal and technical specifications:

- Machines Directive;
- Regulations for the prevention of industrial accidents;
- Workes health and safety improvement at work
- Safety improvement at equipment
- Low voltage directive;
- Electromagnetic compatibility directive.

The manufacturer reserves the material and intellectual property rights inherent to this publication, prohibits its issue and duplication, even partially, without his prior written consent.

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Declares that the BMTECH 4.5 appliance complies with the safety regulations concerning the EC directives for ELECTROMAGNETIC COMPATIBILITY (2014/30 / EU) AND LOW VOLTAGE (2014/35 / EU) and the Machines Directive (2006/42 / EC).

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**FARMALABOR Srl**